

Standard Operating Procedure

Title: Identification of Microorganisms to Genus and Species Level

Note: Avoid picking up too much growth as this will result in a heavy smear, which will not stain effectively and cannot be examined with accuracy.

Note: You may wish to vortex a broth culture prior to creating a smear as this will ensure an even mix of cells in the culture. As the broth is a liquid culture, water is not required to create the smear.

- 3.1.2. Allow the smear to air dry and heat fix the cells to the slide by passing the slide through the hottest part of the flame three times. Allow to cool.
- 3.1.3. Flood the slide with Crystal Violet solution and stand for 30 seconds (Fig. 2A).
- 3.1.4. Rinse off Crystal Violet solution with water and tilt the slide to drain excess water.
- 3.1.5. Flood the slide with 1% Lugol's Iodine and stand for 30 seconds (Fig. 2B).
- 3.1.6. Rinse off Iodine with water and tilt the slide to drain excess water.

- 3.1.7. Flood the slide with Decolourising solution and immediately rinse the slide with running water after 3 seconds (Fig. 2C). Tilt the slide to drain excess water.

*Note: Prolonged contact with decolourising solution will result in over-decolourisation of cells and a false gram reaction. Ensure the contact time is limited to **3 seconds only**.*

- 3.1.8. Flood the slide with the desired counterstaining solution. Leave to act for *up to 2* minutes (Fig. 2D).
- 3.1.9. Rinse with water and blot dry using a paper towel. Do not rub.
- 3.1.10. Examine the smears under Oil Immersion.
- 3.1.11. Record the cellular arrangement (see Appendix 7.2) and Gram stain result of the organism on the Isolate Identification Record form.

3.1.12. Interpretation

- **Gram-positive** Gram-positive organisms stain purple
- **Gram-negative** Gram-negative organisms stain pink/red.

3.1.13. Quality Control

3.1.13.1. Gram-positive and Gram-negative reference cultures are to be used for quality control of the gram staining procedure and reagents.

- **Gram Positive** *Staphylococcus aureus* ATCC 6538
- **Gram Negative** *Escherichia coli* ATCC 8739

3.2. Technical Information

Some Gram-positive bacteria may appear Gram-negative in whole or in part eg. some strains of Gram-variable *Bacillus* and *Clostridium*. Yeasts will appear "Gram-positive" and will be visibly much larger than a prokaryotic bacterial cell.

All Gram-variable and suspected over or under-decolourised organisms are to be subject to a KOH test (see section 3.3).